




01

SIEVES 60mm.Ø x 18mm.H STAINLESS STEEL FRAME

60Ø x 18H.



A	60,4 mm. -0,5
B	60,6 mm. +0,5
C	18 mm. Nominal
D	25 mm. Nominal

- *For small sample sizes
- *Suitable for field work
- *Each sieve is marked with a serial number to ensure traceability
- *Delivered with a Certificate of Compliance according to 2.1 EN 10204, at no extra charge
- *Certificate of Calibration indicating uncertainty of sieve (on request)

WEIGHT OF SIEVES : From 32g to 50g depending on the mesh.

RECOMMENDED VOLUMES OF TEST MATERIAL TO BE SIEVED: Up to 5 cm³.

CONFIGURATION AVAILABLE

01.1.SIEVES WITH METALLIC MESH

STANDARDS: ISO 3310.1 – ASTM E-11-95

01.1.SIEVES WITH METALLIC MESH

ISO 3310.1	ASTM E-11-95	ARTIC.No CISA
20 µ	No 635	0100601
25 µ	No 500	0101001
32 µ	No 450	0101001
36 µ		0101201
38 µ	No 400	0101401
40 µ		0101601
45 µ	No 325	0101801
50 µ		0102001
53 µ	No 270	0102201
56 µ		0102401
63 µ	No 230	0102601
71 µ		0102801
75 µ	No 200	0103001
80 µ		0103201
90 µ	No 170	0103401
100 µ		0103601
106 µ	No 140	0103801
112 µ		0104001
125 µ	No 120	0104201
140 µ		0104401
150 µ	No 100	0104601
160 µ		0104801
180 µ	No 80	0105001
200 µ		0105201
212 µ	No 70	0105401
224 µ		0105601
250 µ	No 60	0105801
280 µ		0106001
300 µ	No 50	0106201
315 µ		0106401
355 µ	No 45	0106601
400 µ		0106801
425 µ	No 40	0107001
450 µ		0107201
500 µ	No 35	0107401
560 µ		0107601
600 µ	No 30	0107801
630 µ		0108001
710 µ	No 25	0108201
800 µ		0108401
850 µ	No 20	0108601
900 µ		0108801

ISO 3310.1	ASTM E-11-95	ARTIC.No CISA
1.00 mm	No 18	0109001
1.12 mm		0109201
1.18 mm	No 16	0109401
1.25 mm		0109601
1.40 mm	No 14	0109801
1.60 mm		0110001
1.70 mm	No 12	0110201
1.80 mm		0110401
2.00 mm	No 10	0110601

ACCESSORIES

DESCRIPTION ARTIC.No

COVER AND PAN

60mm.Ø 01000000

ON REQUEST

● **CERTIFICATE OF CALIBRATION ACCORDING TO:**

ISO 3310

ASTM E-11

AFNOR NF X-11-504

BS 410

● **SIEVES WITH NYLON**

MESH FROM 5µ TO

2.00mm.



msscscientific

Chromatographie-Handel GmbH

Gneisenaustrasse 66/67 · 10961 Berlin · Germany

Fon: 030 / 6270 6087 · Fax: 030 / 6270 6089

E-Mail: info@msscscientific.de · Web: www.msscscientific.de