

## Price List - Laboratory Porcelain, Aluminium Oxide Products, Agate Mortars and Pestles

Part No.	Description	Glaze	Max. Temp.	Conformity	Qty.	Price
<b>Laboratory porcelain</b>						
<b>Crucibles</b>						
LC-J-116	Crucible, tall form, 15mm high, 15mm OD, 1.3ml	G2	1150 °C		each	1,52
LC-J-117	Crucible, tall form, 20mm high, 24mm OD, 5ml	G2	1150 °C		each	1,68
LC-J-118	Crucible, tall form, 26mm high, 31mm OD, 10ml	G2	1150 °C		each	1,76
LC-J-119	Crucible, tall form, 29mm high, 35mm OD, 15ml	G2	1150 °C	ISO 1772, DIN 12 904	each	2,24
LC-J-234	Crucible, tall form, 32mm high, 40mm OD, 25ml	G2	1150 °C	ISO 1772, DIN 12 904	each	2,40
LC-J-120	Crucible, tall form, 37mm high, 43mm OD, 30ml	G2	1150 °C	ISO 1772, DIN 12 904	each	2,48
LC-J-121	Crucible, tall form, 40mm high, 50mm OD, 40ml	G2	1150 °C	ISO 1772, DIN 12 904	each	2,96
LC-J-122	Crucible, tall form, 43mm high, 53mm OD, 50ml	G2	1150 °C		each	3,36
LC-J-123	Crucible, tall form, 54mm high, 65mm OD, 100ml	G2	1150 °C		each	5,12
LC-J-124	Crucible, tall form, 72mm high, 88mm OD, 250ml	G2	1150 °C		each	6,56
LC-J-134	Crucible, wide form, 12mm high, 18mm OD, 2ml	G2	1150 °C		each	1,52
LC-J-517	Crucible, wide form, 16mm high, 26mm OD, 5ml	G2	1150 °C		each	1,76
LC-J-135	Crucible, wide form, 20mm high, 32mm OD, 10ml	G2	1150 °C	ISO 1772, DIN 12 904	each	2,16
LC-J-136	Crucible, wide form, 26mm high, 41mm OD, 20ml	G2	1150 °C	ISO 1772, DIN 12 904	each	2,24
LC-J-137	Crucible, wide form, 31mm high, 50mm OD, 30ml	G2	1150 °C	ISO 1772, DIN 12 904	each	2,96
LC-J-138	Crucible, wide form, 37mm high, 61mm OD, 50ml	G2	1150 °C	ISO 1772, DIN 12 904	each	4,88
LC-J-139	Crucible, wide form, 46mm high, 76mm OD, 100ml	G2	1150 °C		each	5,12
LC-J-140	Crucible, wide form, 52mm high, 84mm OD, 150ml	G2	1150 °C		each	5,76
LC-J-141	Crucible, wide form, 60mm high, 102mm OD, 250ml	G2	1150 °C		each	7,04
LC-J-518	Crucible, wide form, 87mm high, 120mm OD, 500ml	G2	1150 °C		each	10,40
LC-J-143	Gooch crucible, 29mm OD, 15ml, for 15–16mm diameter filter paper	G2	1150 °C		each	2,64
LC-J-144	Gooch crucible, 36mm OD, 30ml, for 20–21mm diameter filter paper	G2	1150 °C		each	4,08
LC-J-145	Gooch crucible, 40mm OD, 50ml, for 23–24mm diameter filter paper	G2	1150 °C		each	5,60
LC-J-146	Gooch crucible, 55mm OD, 140ml, for 30mm diameter filter paper	G2	1150 °C		each	8,80
LC-J-148	Crucible with porous bottom, 36mm OD, 25ml, approx. 1.2µm porosity	G12	autoclavable		each	22,56
LC-J-149	Crucible with porous bottom, 36mm OD, 25ml, approx. 5µm porosity	G12	autoclavable		each	22,56
LC-J-150	Crucible with porous bottom, 36mm OD, 25ml, approx. 15µm porosity	G12	autoclavable		each	22,56
<b>Crucible lids</b>						
LC-J-125	Crucible cover, 26mm OD, for crucible LC-J-117	G3	1150 °C		each	0,96
LC-J-126	Crucible cover, 34mm OD, for crucible LC-J-118, LC-J-134, LC-J-135, LC-J-143	G3	1150 °C		each	1,20
LC-J-127	Crucible cover, 38mm OD, for crucible LC-J-119, LC-J-144, LC-J-147, LC-J-148, LC-J-149, LC-J-150	G3	1150 °C	ISO 1772, DIN 12 904	each	1,04
LC-J-128	Crucible cover, 47mm OD, for crucible LC-J-120, LC-J-136, LC-J-142, LC-J-145	G3	1150 °C	ISO 1772, DIN 12 904	each	1,28
LC-J-129	Crucible cover, 53mm OD, for crucible LC-J-121, LC-J-137	G3	1150 °C	ISO 1772, DIN 12 904	each	1,44
LC-J-529	Crucible cover, 58mm OD, for crucible LC-J-122, LC-J-146	G3	1150 °C		each	1,52
LC-J-130	Crucible cover, 69mm OD, for crucible LC-J-123, LC-J-138	G3	1150 °C		each	1,92
LC-J-131	Crucible cover, 80mm OD, for crucible LC-J-139	G3	1150 °C		each	1,92
LC-J-132	Crucible cover, 93mm OD, for crucible LC-J-124, LC-J-140	G3	1150 °C		each	2,64
LC-J-133	Crucible cover, 108mm OD, for crucible LC-J-141	G3	1150 °C		each	3,12
<b>Evaporating dishes</b>						
LC-J-163	Evaporating dish, standard form, 60mm OD, 35ml	G5	1150 °C		each	3,68
LC-J-164	Evaporating dish, standard form, 75mm OD, 70ml	G5	1150 °C		each	3,68

## Price List - Laboratory Porcelain, Aluminium Oxide Products, Agate Mortars and Pestles

Part No.	Description	Glaze	Max. Temp.	Conformity	Qty.	Price
LC-J-165	Evaporating dish, standard form, 80mm OD, 80ml	G5	1150 °C		each	4,08
LC-J-166	Evaporating dish, standard form, 94mm OD, 120ml	G5	1150 °C		each	5,12
LC-J-167	Evaporating dish, standard form, 100mm OD, 150ml	G5	1150 °C		each	5,84
LC-J-168	Evaporating dish, standard form, 115mm OD, 250ml	G5	1150 °C		each	7,28
LC-J-169	Evaporating dish, standard form, 145mm OD, 385ml	G5	1150 °C		each	12,56
LC-J-170	Evaporating dish, standard form, 162mm OD, 525ml	G5	1150 °C		each	18,48
LC-J-171	Evaporating dish, standard form, 185mm OD, 765ml	G5	1150 °C		each	28,16
LC-J-172	Evaporating dish, standard form, 215mm OD, 1285ml	G5	1150 °C		each	33,76
LC-J-173	Evaporating dish, standard form, 265mm OD, 2100ml	G5	1150 °C		each	53,60
LC-J-174	Evaporating dish, standard form, 305mm OD, 3250ml	G5	1150 °C		each	61,52
LC-J-613	Evaporating dish, standard form, 200mm OD, 1700ml	G5	1150 °C		each	41,20
LC-J-175	Evaporating dish, low form, 77mm OD, 100ml	G12	autoclavable		each	3,68
LC-J-176	Evaporating dish, low form, 65mm OD, 35ml	G6	1150 °C		each	2,72
LC-J-177	Evaporating dish, low form, 75mm OD, 50ml	G6	1150 °C		each	2,88
LC-J-178	Evaporating dish, low form, 85mm OD, 80ml	G6	1150 °C		each	3,44
LC-J-179	Evaporating dish, low form, 100mm OD, 100ml	G6	1150 °C		each	4,40
LC-J-180	Evaporating dish, low form, 120mm OD, 200ml	G6	1150 °C		each	5,28
LC-J-681	Evaporating dish, low form, 150mm OD, 300ml	G6	1150 °C		each	10,88
LC-J-682	Evaporating dish, low form, 175mm OD, 500ml	G6	1150 °C		each	15,44
LC-J-683	Evaporating dish, low form, 200mm OD, 1000ml	G6	1150 °C		each	26,00
LC-J-684	Evaporating dish, low form, 250mm OD, 1500ml	G6	1150 °C		each	51,36
LC-J-685	Evaporating dish, low form, 300mm OD, 2000ml	G6	1150 °C		each	94,72
LC-J-105	Evaporating dish with handle (casserole), 40mm OD, 25ml	G4	1150 °C		each	3,92
LC-J-106	Evaporating dish with handle (casserole), 70mm OD, 65ml	G4	1150 °C		each	5,28
LC-J-107	Evaporating dish with handle (casserole), 85mm OD, 140ml	G4	1150 °C		each	7,04
LC-J-108	Evaporating dish with handle (casserole), 100mm OD, 200ml	G4	1150 °C		each	8,80
LC-J-109	Evaporating dish with handle (casserole), 125mm OD, 500ml	G4	1150 °C		each	15,68
LC-J-110	Evaporating dish with handle (casserole), 150mm OD, 760ml	G4	1150 °C		each	21,60
LC-J-111	Evaporating dish with handle (casserole), 175mm OD, 1250ml	G4	1150 °C		each	26,48
<b>Mortars and pestles</b>						
LC-J-701	Mortar & pestle, 80mm OD, 60ml	G7	autoclavable		each	5,76
LC-J-702	Mortar & pestle, 90mm OD, 100ml	G7	autoclavable		each	7,44
LC-J-703	Mortar & pestle, 105mm OD, 150ml	G7	autoclavable		each	10,24
LC-J-704	Mortar & pestle, 125mm OD, 275ml	G7	autoclavable		each	14,80
LC-J-705	Mortar & pestle, 155mm OD, 550ml	G7	autoclavable		each	23,92
LC-J-706	Mortar & pestle, 180mm OD, 900ml	G7	autoclavable		each	31,52
LC-J-707	Mortar & pestle, 200mm OD, 1500ml	G7	autoclavable		each	53,76
LC-J-708	Mortar & pestle, 230mm OD, 2000ml	G7	autoclavable		each	63,04
LC-J-709	Mortar & pestle, 285mm OD, 4200ml	G7	autoclavable		each	178,64
LC-J-710	Mortar & pestle, 330mm OD, 7000ml	G7	autoclavable		each	273,44
LC-J-198	Mortar & pestle, deep form, 70mm OD, 50ml	G7	autoclavable		each	6,08
LC-J-199	Mortar & pestle, deep form, 80mm OD, 70ml	G7	autoclavable		each	6,80
LC-J-200	Mortar & pestle, deep form, 90mm OD, 150ml	G7	autoclavable		each	9,04

## Price List - Laboratory Porcelain, Aluminium Oxide Products, Agate Mortars and Pestles

Part No.	Description	Glaze	Max. Temp.	Conformity	Qty.	Price
LC-J-201	Mortar & pestle, deep form, 115mm OD, 275ml	G7	autoclavable		each	11,44
LC-J-202	Mortar & pestle, deep form, 130mm OD, 400ml	G7	autoclavable		each	16,72
LC-J-203	Mortar & pestle, deep form, 165mm OD, 750ml	G7	autoclavable		each	24,00
LC-J-204	Mortar & pestle, deep form, 210mm OD, 1900ml	G7	autoclavable		each	52,32
LC-J-205	Mortar & pestle, deep form, 285mm OD, 4000ml	G7	autoclavable		each	125,28
LC-J-214	Mortar & pestle, deep form, 256mm OD, 2750ml	G7	autoclavable		each	89,04
LC-J-252	Pestle, length 93mm, for mortar LC-J-701, LC-J-702	G8	autoclavable		each	2,72
LC-J-206	Pestle, length 114mm, for mortar LC-J-703, LC-J-198	G8	autoclavable		each	3,44
LC-J-207	Pestle, length 130mm, for mortar LC-J-704, LC-J-199	G8	autoclavable		each	4,32
LC-J-208	Pestle, length 157mm, for mortar LC-J-705, LC-J-200	G8	autoclavable		each	5,76
LC-J-209	Pestle, length 180mm, for mortar LC-J-706, LC-J-201	G8	autoclavable		each	7,04
LC-J-210	Pestle, length 194mm, for mortar LC-J-707, LC-J-202	G8	autoclavable		each	8,56
LC-J-211	Pestle, length 222mm, for mortar LC-J-708, LC-J-203	G8	autoclavable		each	13,36
LC-J-212	Pestle, length 262mm, for mortar LC-J-709, LC-J-204	G8	autoclavable		each	18,96
LC-J-213	Pestle, length 280mm, for mortar LC-J-710, LC-J-205, LC-J-214	G8	autoclavable		each	25,52
<b>Funnels</b>						
LC-J-183	Buchner funnel, 20mm OD, 3ml, for 13mm diameter filter paper	G3	autoclavable		each	3,84
LC-J-184	Buchner funnel, 50mm OD, 35ml, for 35–40mm diameter filter paper	G3	autoclavable		each	6,56
LC-J-185	Buchner funnel, 65mm OD, 100ml, for 50–55mm diameter filter paper	G3	autoclavable		each	9,68
LC-J-186	Buchner funnel, 86mm OD, 200ml, for 70mm diameter filter paper	G3	autoclavable		each	10,48
LC-J-187	Buchner funnel, 105mm OD, 350ml, for 90mm diameter filter paper	G3	autoclavable		each	13,36
LC-J-188	Buchner funnel, 127mm OD, 600ml, for 110mm diameter filter paper	G3	autoclavable		each	21,04
LC-J-189	Buchner funnel, 143mm OD, 800ml, for 110–124mm diameter filter paper	G3	autoclavable		each	24,56
LC-J-190	Buchner funnel, 172mm OD, 1150ml, for 142–150mm diameter filter paper	G3	autoclavable		each	50,72
LC-J-191	Buchner funnel, 210mm OD, 2000ml, for 185mm diameter filter paper	G3	autoclavable		each	78,08
LC-J-192	Buchner funnel, 270mm OD, 4800ml, for 240mm diameter filter paper	G3	autoclavable		each	223,52
LC-J-780	Hirsch funnel, 30mm OD, 10ml, for 10mm diameter filter paper	G12	autoclavable		each	4,32
LC-J-193	Hirsch funnel, 40mm OD, 12ml, for 13–15mm diameter filter paper	G12	autoclavable		each	5,12
LC-J-195	Hirsch funnel, 50mm OD, 25ml, for 25–30mm diameter filter paper	G12	autoclavable		each	5,84
LC-J-196	Hirsch funnel, 75mm OD, 50ml, for 40–47mm diameter filter paper	G12	autoclavable		each	9,20
LC-J-197	Hirsch funnel, 100mm OD, 140ml, for 55mm diameter filter paper	G12	autoclavable		each	10,24
<b>Combustion boats</b>						
LC-J-801	Combustion boat, length 80mm, 1,5ml	G1	1150 °C		each	2,40
LC-J-802	Combustion boat, length 86mm, 2ml	G1	1150 °C		each	2,48
LC-J-803	Combustion boat, length 90mm, 7ml	G1	1150 °C		each	2,64
LC-J-804	Combustion boat, length 120mm, 12ml	G1	1150 °C		each	3,20
LC-J-801-U	Combustion boat, unglazed, length 80mm, 1,5ml	G12	1150 °C		each	1,44
LC-J-802-U	Combustion boat, unglazed, length 86mm, 2ml	G12	1150 °C		each	1,68
LC-J-803-U	Combustion boat, unglazed, length 90mm, 7ml	G12	1150 °C		each	1,68
LC-J-804-U	Combustion boat, unglazed, length 120mm, 12ml	G12	1150 °C		each	2,24
<b>Dessicator plates</b>						
LC-J-901	Dessicator plate, without feet, 230mm OD, 8 x 30mm Ø holes	G9	autoclavable		each	15,68
LC-J-931	Dessicator plate, without feet, 95mm OD, 3 x 30mm Ø holes	G9	autoclavable		each	10,16

## Price List - Laboratory Porcelain, Aluminium Oxide Products, Agate Mortars and Pestles

Part No.	Description	Glaze	Max. Temp.	Conformity	Qty.	Price
LC-J-932	Dessicator plate, without feet, 115mm OD, 3 x 30mm Ø holes	G9	autoclavable		each	12,32
LC-J-933	Dessicator plate, without feet, 140mm OD, 6 x 30mm Ø holes	G9	autoclavable		each	14,08
LC-J-934	Dessicator plate, without feet, 165mm OD, 7 x 30mm Ø holes	G9	autoclavable		each	15,92
LC-J-935	Dessicator plate, without feet, 190mm OD, 8 x 30mm Ø holes	G9	autoclavable		each	20,32
LC-J-936	Dessicator plate, without feet, 230mm OD, 10 x 30mm Ø holes	G9	autoclavable		each	21,84
LC-J-937	Dessicator plate, without feet, 290mm OD, 13 x 30mm Ø holes	G9	autoclavable		each	30,00
LC-J-153	Dessicator plate, three small feet, 95mm AD, 3 x 30mm Ø holes	G9	1150 °C		each	12,08
LC-J-154	Dessicator plate, three small feet, 115mm AD, 4 x 30mm Ø holes	G9	1150 °C		each	16,00
LC-J-155	Dessicator plate, three small feet, 142mm AD, 5 x 30mm Ø holes	G9	1150 °C		each	16,00
LC-J-156	Dessicator plate, three small feet, 145mm AD, 7 x 30mm Ø holes	G9	1150 °C		each	18,40
LC-J-157	Dessicator plate, three small feet, 190mm AD, 7 x 30mm Ø holes	G9	1150 °C		each	21,60
LC-J-158	Dessicator plate, three small feet, 230mm AD, 8 x 30mm Ø holes	G9	1150 °C		each	23,44
LC-J-159	Dessicator plate, without feet, numerous 5mm Ø holes, centre hole 23mm Ø, 140mm OD	G9	1150 °C		each	16,72
LC-J-160	Dessicator plate, without feet, numerous 5mm Ø holes, centre hole 23mm Ø, 190mm OD	G9	1150 °C		each	18,64
LC-J-161	Dessicator plate, without feet, numerous 5mm Ø holes, centre hole 23mm Ø, 230mm OD	G9	1150 °C		each	21,84
LC-J-1150	Dessicator plate, without feet, numerous 5mm Ø holes, centre hole 23mm Ø, 290mm OD	G9	1150 °C		each	28,56
LC-J-162	Dessicator plate, without feet, numerous 2mm Ø holes, centre hole 23mm Ø, 230mm OD	G9	1150 °C		each	21,60
<b>Spatulas</b>						
LC-J-225	Spatula, length 205mm	G10	autoclavable		each	10,96
LC-J-226	Spatula with spoon, length 100mm, 0,1ml	G10	autoclavable		each	3,68
LC-J-227	Spatula with spoon, length 125mm, 0,3ml	G10	autoclavable		each	4,64
LC-J-228	Spatula with spoon, length 165mm, 0,5ml	G10	autoclavable		each	5,36
LC-J-229	Spatula with spoon, length 198mm, 1,2ml	G10	autoclavable		each	8,48
LC-J-230	Spatula with spoon, length 260mm, 5,5ml	G10	autoclavable		each	10,48
<b>Spotting-Plates/Streak plate</b>						
LC-J-820	Spotting-/Color-Plate, 94 x 62mm, 6 depressions	G1	autoclavable		each	4,80
LC-J-221	Spotting-/Color-Plate, 115 x 88mm, 12 depressions	G1	autoclavable		each	7,68
LC-J-223	Spotting-/Color-Plate, 117 x 137mm, 30 depressions	G1	autoclavable		each	22,72
LC-J-224	Streak plate, 65 x 50 x 3mm	G12	autoclavable		each	2,96
<b>Clay pipe triangles</b>						
LC-J-884	Clay pipe triangle, side length 50mm	G12	autoclavable		each	1,68
LC-J-885	Clay pipe triangle, side length 60mm	G12	autoclavable		each	1,76
LC-J-886	Clay pipe triangle, side length 75mm	G12	autoclavable		each	2,00
<b>Beehive Shelves</b>						
LC-J-651	Beehive Shelf, 60mm OD	G11	autoclavable		each	2,16
LC-J-652	Beehive Shelf, 75mm OD	G11	autoclavable		each	2,40
LC-J-653	Beehive Shelf, 85mm OD	G11	autoclavable		each	2,96
LC-J-654	Beehive Shelf, 100mm OD	G11	autoclavable		each	4,32
<b>Porous pots</b>						
LC-J-874	Porous pot, 50mm OD, 150ml				each	2,88
LC-J-875	Porous pot, 60mm OD, 150ml				each	3,12
LC-J-876	Porous pot, 75mm OD, 150ml				each	4,88

## Price List - Laboratory Porcelain, Aluminium Oxide Products, Agate Mortars and Pestles

Part No.	Description	Glaze	Max. Temp.	Conformity	Qty.	Price
<b>Aluminium oxide products</b>						
<b>Aluminium oxide crucibles</b>						
LC-HACC.0.2	Crucible, cylindrical, 6mm high, 10mm OD, 0.2ml		1750 °C		each	8,13
LC-HACC-1	Crucible, cylindrical, 12mm high, 17mm OD, 1ml		1750 °C		each	9,85
LC-HACC-3	Crucible, cylindrical, 22mm high, 17mm OD, 3ml		1750 °C		each	11,58
LC-HACC-5	Crucible, cylindrical, 47mm high, 17mm OD, 5ml		1750 °C		each	14,15
LC-HACC-5-A	Crucible, cylindrical, 17mm high, 47mm OD, 5ml		1750 °C		each	14,15
LC-HACC-10	Crucible, cylindrical, 29mm high, 28mm OD, 10ml		1750 °C		each	15,38
LC-HACC-15	Crucible, cylindrical, 42mm high, 28mm OD, 15ml		1750 °C		each	18,40
LC-HACC-15-A	Crucible, cylindrical, 35mm high, 32mm OD, 15ml		1750 °C		each	18,40
LC-HACC-20	Crucible, cylindrical, 40mm high, 32mm OD, 20ml		1750 °C		each	21,43
LC-HACC-20-A	Crucible, cylindrical, 50mm high, 28mm OD, 20ml		1750 °C		each	21,43
LC-HACC-30	Crucible, cylindrical, 45mm high, 35mm OD, 30ml		1750 °C		each	25,30
LC-HACC-40	Crucible, cylindrical, 52mm high, 40mm OD, 40ml		1750 °C		each	28,73
LC-HACC-50	Crucible, cylindrical, 56mm high, 41mm OD, 50ml		1750 °C		each	29,98
LC-HACC-60	Crucible, cylindrical, 40mm high, 55mm OD, 60ml		1750 °C		each	32,90
LC-HACC-60-A	Crucible, cylindrical, 74mm high, 41mm OD, 60ml		1750 °C		each	32,90
LC-HACC-80	Crucible, cylindrical, 65mm high, 47mm OD, 80ml		1750 °C		each	37,25
LC-HACC-80-A	Crucible, cylindrical, 93mm high, 40mm OD, 80ml		1750 °C		each	37,25
LC-HACC-100	Crucible, cylindrical, 65mm high, 52mm OD, 100ml		1750 °C		each	41,45
LC-HACC-150	Crucible, cylindrical, 89mm high, 54mm OD, 150ml		1750 °C		each	54,33
LC-HACC-150-A	Crucible, cylindrical, 55mm high, 70mm OD, 150ml		1750 °C		each	54,33
LC-HACC-170	Crucible, cylindrical, 87mm high, 57mm OD, 170ml		1750 °C		each	62,03
LC-HACC-225	Crucible, cylindrical, 83mm high, 64mm OD, 225ml		1750 °C		each	71,50
LC-HACC-250	Crucible, cylindrical, 68mm high, 80mm OD, 250ml		1750 °C		each	80,00
LC-HACC-275	Crucible, cylindrical, 75mm high, 80mm OD, 275ml		1750 °C		each	96,58
LC-HACC-410	Crucible, cylindrical, 82mm high, 90mm OD, 410ml		1750 °C		each	132,03
LC-HACC-450	Crucible, cylindrical, 105mm high, 84mm OD, 450ml		1750 °C		each	155,68
LC-HACC-750	Crucible, cylindrical, 215mm high, 85mm OD, 750ml		1750 °C		each	269,48
LC-HACC-1000	Crucible, cylindrical, 168mm high, 100mm OD, 1000ml		1750 °C		each	385,03
LC-HACC-1500	Crucible, cylindrical, 193mm high, 105mm OD, 1500ml		1750 °C		each	547,53
LC-HACC-2500	Crucible, cylindrical, 247mm high, 127mm OD, 2500ml		1750 °C		each	1.026,80
LC-HACHF-5	Crucible, tall form, 29mm high, 23mm OD, 5ml		1750 °C		each	14,15
LC-HACHF-10	Crucible, tall form, 23mm high, 30mm OD, 10ml		1750 °C		each	15,38
LC-HACHF-20	Crucible, tall form, 36mm high, 39mm OD, 20ml		1750 °C		each	21,43
LC-HACHF-50	Crucible, tall form, 47mm high, 52mm OD, 50ml		1750 °C		each	29,98
LC-HACHF-100	Crucible, tall form, 65mm high, 62mm OD, 100ml		1750 °C		each	45,33
LC-HACHF-250	Crucible, tall form, 90mm high, 82mm OD, 250ml		1750 °C		each	108,20
LC-HACLF-10	Crucible, wide form, 25mm high, 39mm OD, 10ml		1750 °C		each	15,38
LC-HACLF-20	Crucible, wide form, 35mm high, 50mm OD, 20ml		1750 °C		each	21,43
LC-HACLF-25	Crucible, wide form, 37mm high, 43mm OD, 25ml		1750 °C		each	24,38
LC-HACLF-30	Crucible, wide form, 40mm high, 56mm OD, 30ml		1750 °C		each	25,30
LC-HACLF-50	Crucible, wide form, 44mm high, 61mm OD, 50ml		1750 °C		each	29,98

## Price List - Laboratory Porcelain, Aluminium Oxide Products, Agate Mortars and Pestles

Part No.	Description	Glaze	Max. Temp.	Conformity	Qty.	Price
LC-HACLF-80	Crucible, wide form, 56mm high, 64mm OD, 80ml		1750 °C		each	37,25
LC-HACLF-175	Crucible, wide form, 50mm high, 82mm OD, 175ml		1750 °C		each	64,48
LC-HACLF-225	Crucible, wide form, 58mm high, 87mm OD, 225ml		1750 °C		each	82,58
LC-HACCO-2	Crucible, conical, 18mm high, 18mm OD, 2ml		1750 °C		each	10,70
LC-HACCO-5	Crucible, conical, 35mm high, 27mm OD, 5ml		1750 °C		each	12,88
LC-HACCO-20	Crucible, conical, 42mm high, 33mm OD, 20ml		1750 °C		each	20,98
LC-HACCO-30	Crucible, conical, 43mm high, 48mm OD, 30ml		1750 °C		each	28,73
LC-HACCO-50	Crucible, conical, 53mm high, 46mm OD, 50ml		1750 °C		each	33,85
LC-HACCO-100	Crucible, conical, 86mm high, 64mm OD, 100ml		1750 °C		each	62,48
LC-HACCO-150	Crucible, conical, 90mm high, 63mm OD, 150ml		1750 °C		each	50,90
LC-HACCO-170	Crucible, conical, 61mm high, 80mm OD, 170ml		1750 °C		each	73,18
LC-HACCO-200	Crucible, conical, 94mm high, 73mm OD, 200ml		1750 °C		each	98,40
LC-HACCO-250	Crucible, conical, 100mm high, 82mm OD, 250ml		1750 °C		each	115,63
LC-HACCO-300	Crucible, conical, 87mm high, 91mm OD, 300ml		1750 °C		each	115,90
LC-HACCO-500	Crucible, conical, 128mm high, 88mm OD, 500ml		1750 °C		each	172,75
LC-HACCO-750	Crucible, conical, 125mm high, 125mm OD, 750ml		1750 °C		each	334,38
<b>Aluminium oxide crucible lids</b>						
LC-HACD-10	Crucible lid, 10mm OD		1750 °C		each	6,53
LC-HACD-15	Crucible lid, 15mm OD		1750 °C		each	8,13
LC-HACD-20	Crucible lid, 20mm OD		1750 °C		each	10,55
LC-HACD-25	Crucible lid, 25mm OD		1750 °C		each	10,70
LC-HACD-28	Crucible lid, 28mm OD		1750 °C		each	12,43
LC-HACD-32	Crucible lid, 32mm OD		1750 °C		each	12,43
LC-HACD-38	Crucible lid, 38mm OD		1750 °C		each	14,15
LC-HACD-48	Crucible lid, 48mm OD		1750 °C		each	16,68
LC-HACD-58	Crucible lid, 58mm OD		1750 °C		each	19,25
LC-HACD-75	Crucible lid, 75mm OD		1750 °C		each	24,38
LC-HACD-90	Crucible lid, 90mm OD		1750 °C		each	32,90
LC-HACD-100	Crucible lid, 100mm OD		1750 °C		each	35,55
LC-HACD-115	Crucible lid, 115mm OD		1750 °C		each	50,90
LC-HACD-125	Crucible lid, 125mm OD		1750 °C		each	61,65
LC-HACD-150	Crucible lid, 150mm OD		1750 °C		each	76,60
<b>Aluminium oxide evaporating dishes</b>						
LC-HART-5	Rectangular dish, 50 x 26 x 10mm, 5ml		1750 °C		each	22,65
LC-HART-10	Rectangular dish, 50 x 25 x 19mm, 10ml		1750 °C		each	25,30
LC-HART-10A	Rectangular dish, 98 x 25 x 10mm, 10ml		1750 °C		each	28,73
LC-HART-10B	Rectangular dish, 50 x 50 x 15mm, 10ml		1750 °C		each	32,90
LC-HART-15	Rectangular dish, 85 x 34 x 12mm, 15ml		1750 °C		each	33,85
LC-HART-25	Rectangular dish, 99 x 25 x 20mm, 25ml		1750 °C		each	37,25
LC-HART-25A	Rectangular dish, 77 x 51 x 12mm, 25ml		1750 °C		each	37,25
LC-HART-30	Rectangular dish, 50 x 50 x 27mm, 30ml		1750 °C		each	39,75
LC-HART-80	Rectangular dish, 124 x 40 x 26mm, 80ml		1750 °C		each	53,03
LC-HART-80A	Rectangular dish, 80 x 55 x 30mm, 80ml		1750 °C		each	53,03

## Price List - Laboratory Porcelain, Aluminium Oxide Products, Agate Mortars and Pestles

Part No.	Description	Glaze	Max. Temp.	Conformity	Qty.	Price
LC-HART-100	Rectangular dish, 101 x 52 x 32mm, 100ml		1750 °C		each	69,78
LC-HART-100A	Rectangular dish, 125 x 25 x 40mm, 100ml		1750 °C		each	69,78
LC-HART-130	Rectangular dish, 150 x 67 x 19mm, 130ml		1750 °C		each	105,70
LC-HART-130A	Rectangular dish, 130 x 60 x 25mm, 130ml		1750 °C		each	116,80
LC-HART-140	Rectangular dish, 93 x 52 x 42mm, 140ml		1750 °C		each	116,80
LC-HART-160	Rectangular dish, 120 x 120 x 20mm, 160ml		1750 °C		each	152,98
LC-HART-225	Rectangular dish, 127 x 51 x 40mm, 225ml		1750 °C		each	182,68
LC-HART-240	Rectangular dish, 132 x 52 x 42mm, 240ml		1750 °C		each	208,25
LC-HART-250	Rectangular dish, 100 x 100 x 40mm, 250ml		1750 °C		each	216,75
LC-HART-340	Rectangular dish, 200 x 90 x 25mm, 340ml		1750 °C		each	253,73
LC-HART-1000	Rectangular dish, 250 x 160 x 40mm, 1000ml		1750 °C		each	415,08
LC-HART-2000	Rectangular dish, 250 x 250 x 50mm, 2000ml		1750 °C		each	641,80
LC-HACDS-10	Circular dish, 10mm high, 40mm OD, 10ml		1750 °C		each	15,38
LC-HACDS-20	Circular dish, 12mm high, 50mm OD, 20ml		1750 °C		each	21,43
LC-HACDS-40	Circular dish, 20mm high, 60mm OD, 40ml		1750 °C		each	29,98
LC-HACDS-100	Circular dish, 32mm high, 75mm OD, 100ml		1750 °C		each	50,08
LC-HACDS-250	Circular dish, 43mm high, 101mm OD, 250ml		1750 °C		each	80,00
LC-HACDS-450	Circular dish, 62mm high, 111mm OD, 450ml		1750 °C		each	141,28
LC-HACDS-750	Circular dish, 40mm high, 165mm OD, 750ml		1750 °C		each	191,23
<b>Aluminium oxide combustion boats</b>						
LC-HACBT-3	Combustion boat with handle and hole, 91 x 12 x 11mm, 3ml		1750 °C		each	16,68
LC-HACBT-7	Combustion boat with handle and hole, 86 x 15 x 13mm, 7ml		1750 °C		each	21,43
LC-HACBT-15	Combustion boat with handle and hole, 110 x 20 x 18mm, 15ml		1750 °C		each	22,65
LC-HACBT-30	Combustion boat with handle and hole, 117 x 30 x 19mm, 30ml		1750 °C		each	32,90
LC-HACBT-30A	Combustion boat with handle and hole, 86 x 27 x 27mm, 30ml		1750 °C		each	32,90
LC-HACBT-40	Combustion boat with handle and hole, 136 x 27 x 21mm, 40ml		1750 °C		each	39,75
LC-HACBT-500	Combustion boat with handle and hole, 255 x 90 x 40mm, 500ml		1750 °C		each	359,40
<b>Aluminium oxide mortars and pestles</b>						
LC-HAM-1.5	Mortar & pestle, 18mm high, 35mm OD, 1.5ml		1750 °C		each	77,60
LC-HAM-7	Mortar & pestle, 26mm high, 50mm OD, 7ml		1750 °C		each	85,38
LC-HAM-15	Mortar & pestle, 31mm high, 65mm OD, 15ml		1750 °C		each	108,65
LC-HAM-30	Mortar & pestle, 37mm high, 75mm OD, 30ml		1750 °C		each	139,70
LC-HAM-80	Mortar & pestle, 47mm high, 100mm OD, 80ml		1750 °C		each	232,83
<b>Agate mortars and pestles</b>						
LC-C-100	Mortar & pestle, upper rim 25mm OD, 20mm ID. Overall height 19mm, interior depth 6mm, 1,5ml				each	52,83
LC-C-101	Mortar & pestle, upper rim 38mm OD, 29mm ID. Overall height 25mm, interior depth 9mm, 8ml				each	56,80
LC-C-102	Mortar & pestle, upper rim 50mm OD, 38mm ID. Overall height 28mm, interior depth 11mm, 10ml				each	61,15
LC-C-103	Mortar & pestle, upper rim 64mm OD, 50mm ID. Overall height 31mm, interior depth 13mm, 12ml				each	72,98
LC-C-104	Mortar & pestle, upper rim 76mm OD, 63mm ID. Overall height 35mm, interior depth 15mm, 22ml				each	105,53
LC-C-105	Mortar & pestle, upper rim 89mm OD, 76mm ID. Overall height 38mm, interior depth 18mm, 50ml				each	178,70
LC-C-106	Mortar & pestle, upper rim 100mm OD, 88mm ID. Overall height 41mm, interior depth 23mm, 70ml				each	211,05
LC-C-107	Mortar & pestle, upper rim 115mm OD, 105mm ID. Overall height 46mm, interior depth 27mm, 120ml				each	327,43
LC-C-108	Mortar & pestle, upper rim 127mm OD, 114mm ID. Overall height 51mm, interior depth 29mm, 135ml				each	406,30

### Price List - Laboratory Porcelain, Aluminium Oxide Products, Agate Mortars and Pestles

Part No.	Description	Glaze	Max. Temp.	Conformity	Qty.	Price
LC-C-109	Mortar & pestle, upper rim 139mm OD, 122mm ID. Overall height 56mm, interior depth 30mm, 170ml				each	669,23
LC-C-110	Mortar & pestle, upper rim 150mm OD, 139mm ID. Overall height 65mm, interior depth 39mm, 290ml				each	834,33

#### Glaze

- G1: complete except underside of base
- G2: complete except underside of base and rim
- G3: complete except rim
- G4: complete except rim and top of handle
- G5: inside completely, outside partially
- G6: inside and outside completely, underside of base partially
- G7: outside of mortar, pestle except grinding surface
- G8: except grinding surface
- G9: one side
- G10: complete
- G11: complete except base and rim
- G12: none