

**Heraeus Noblelight Line Sources**

Type/Description	Part No.	Output	Ignition Voltage	Operating Voltage	Current	Price
<b>High-Pressure Hg Line Sources</b>						
HPK 125W	80017490	125 W	180 V	125 V	1.15 A	1.076,00
<b>Medium-Pressure Hg Line Source</b>						
St 46	80001954	33 W	330 V	55 V ( $\pm 5$ V)	0.6 A	443,00
St 49/1	45000124	35 W	230 V	40 V ( $\pm 5$ V)	1.0 A	373,00
St 75	56001171	80 W	800 V	88 V ( $\pm 7$ V)	0.9 A	527,00
<b>Low-Pressure Hg Line Sources</b>						
93123	80016007	12 W			0.9 A	284,00
HG2	80017485	2.7 W	5 - 7 VDC	3 - 4.5 V	0.15 A	479,00
NK 2/1	56001012	2 W	650 VAC or 900 VDC	220 V ( $\pm 20$ V)	0.01 A	284,00
NK 4/1	56001013	3.5 W	1000 VAC or 900 VDC	200 V ( $\pm 40$ V)	0.02 A	331,00
NK 4/4	56001636	3.5 W	1000 VAC or 900 VDC	220 V ( $\pm 20$ V)	0.02 A	303,00
NK 6/12	56001781	6 W	1000 VAC	320 V	0.025 A	488,00
NK 7/8	80001987	5 W	900 VAC	230 V ( $\pm 30$ V)	0.02 A	431,00
TUV 6WE germicidal lamp	80017693	0.12 W		220 V	0.027 A	80,00
<b>Low-Pressure Na Line Sources</b>						
93122	80048310	14 W	200 V	50 V	0.933 A	334,00
<b>Black Light Blue Lamps</b>						
TW6W	80017695	0.05 W		220 V	0.027 A	109,00
<b>Spectral Calibration Line Sources (Compact Systems)</b>						
C-100S line source compact system	80017392					1.299,00
C-150 lamp holder w/cable and plug	80017393					321,00
C-155-Xe xenon replacement lamp for C100S	80017405					132,00
C-155-Ne neon replacement lamp for C100S	80017404					132,00
C-155-Ar argon replacement lamp for C100S	80017400					132,00
C-155-Kr krypton replacement lamp for C100S	80017403					132,00
C-155-Hg/I mercury/iodine replacement lamp for C100S	80017402					132,00
C-155-Hg mercury replacement lamp for C100S	80017401					132,00