

All prices are in Euros per piece and are subject to change without further notice.  
All brand names and registered trademarks are the property of their respective owners.

## Heraeus Noblelight PID-Lamps

Type/Description	Part No.	Fill Gas	Photon Energy (eV)	Operating Current (mA)	Ignition Time (sec.)	Dimensions (mm)	Minimum Purchase	Price/Each
<b>Photoionisation Lamps, DC Excited</b>								
PAL 118	80017617	argon	11,8	≤0,1	≤2	53,5 x 35,0	5	1.151,00
PAS 118	80017619	argon	11,8	≤0,1	≤2	53,5 x 19,6	5	1.151,00
PKA 106	80017630	krypton	10,6	≤2	≤2	53,5 x 19,6	3	478,00
PKL 106	80017631	krypton	10,6	≤2	≤2	53,5 x 35,0	3	513,00
PKS 106	80017644	krypton	10,6	≤2	≤2	53,5 x 19,6	1	478,00
PXL 096	80017648	xenon	9,6	≤2	≤2	53,5 x 35,0	3	513,00
PXS 096	80017652	xenon	9,6	≤2	≤2	53,5 x 19,6	3	478,00
<b>Photoionisation Lamps, RF Excited</b>								
PAR 118	80017618	argon	11,8	200	1-10	53,0 x 12,7	5	800,00
PHR 067	80017627	hydrogen	6,7	200	1-10	53,0 x 12,7	3	248,00
PHR 067-6	80085951	hydrogen	6,7	26	1-10	30,0 x 6,0	3	248,00
PKR 100	80081014	krypton	10	200	1-10	53,0 x 12,7	3	285,00
PKR 100-6-14	80140698	krypton	10	26	1-10	14,0 x 6,0	3	285,00
PKR 106	80014420	krypton	10,6	200	1-10	53,0 x 12,7	1	248,00
PKR 106-6	80121540	krypton	10,6	26	1-10	30,0 x 6,0	3	248,00
PKR 106-6-14	80092804	krypton	10,6	26	1-10	14,0 x 6,0	3	248,00
PXR 096	80017650	xenon	9,6	200	1-10	53,0 x 12,7	3	248,00

Please note that the PAL, PAR and PAS lamp types should not be stored unpacked for more than 4 weeks (packed approx. 6 months). Their lithium fluoride windows are very hygroscopic.