

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information						Lamp Information											
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price				
ABI	refer to Applied Biosystems																
Abimed	refer to Gilson																
Agilent Technologies	1050 (79853C) VWD	79883-60002									HP-404	J, L	436,00				
Agilent Technologies	1050 (79854A)	79883-60002									HP-404	J, L	436,00				
Agilent Technologies	1050 (G1306A) DAD	79883-60002									HP-404	J, L	436,00				
Agilent Technologies	1090 series I DAD (79880A)	79883-60002									HP-404	J, L	436,00				
Agilent Technologies	1090 series II DAD (79883A)	79883-60002									HP-404	J, L	436,00				
Agilent Technologies	1100 VWD (G1314A)	G1314-60100									DX 224/05 J	80034059	J, L	493,00			
Agilent Technologies	1100 DAD (G1315A)	2140-0590 / 2140-0813 / 5181-1530									LC-D2-AGI-G1315	J, L	321,74	HP-404	J, L	436,00	
Agilent Technologies	1100 MWD (G1365A)	2140-0590 / 2140-0813 / 5181-1530									LC-D2-AGI-G1315	J, L	321,74	HP-404	J, L	436,00	
Agilent Technologies	1200 VWD (G1314B)	G1314-60100									DX 224/05 J	80034059	J, L	493,00	HP-404	J, L	436,00
Agilent Technologies	1200 VWD-SL (G1314C)	G1314-60100									DX 224/05 J	80034059	J, L	493,00	HP-404	J, L	436,00
Agilent Technologies	1200 DAD (G1315B)	2140-0590 / 2140-0813 / 5181-1530									LC-D2-AGI-G1315	J, L	321,74	HP-404	J, L	436,00	
Agilent Technologies	1200 MWD (G1365B)	2140-0590 / 2140-0813 / 5181-1530									LC-D2-AGI-G1315	J, L	321,74	HP-404	J, L	436,00	
Agilent Technologies	1260 DAD										LC-D2-AGI-G1315	J, L, nR	321,74	HP-404	J, L, nR	436,00	
Agilent Technologies	1260 VWD										LC-D2-AGI-G1315	J, L, nR	321,74	HP-404	J, L, nR	436,00	
Agilent Technologies	1290 DAD	5190-0917												LC-D22K-AGI-03	J, L, nR	450,00	
Agilent Technologies	8452A	08452-60104									DS 204/05 RTJ	45000030	J, T	467,00	HP-405	J, L	436,00
Agilent Technologies	8453A	2140-0605												HP-404	J, L	436,00	
Agilent Technologies	CE G 1600 Capillary Electrophoresis System													HP-404	J, L	436,00	
Agilent Technologies	HP 3D CE (G1600A CE)													HP-404	J, L	436,00	
Alpkern	510		J53	80078578	J, T	441,00											
Alpkern	5010		J53	80078578	J, T	441,00											
Altex	100		R07	80013769		416,00	SD 3651-03		80056123	452,00							
Altex	153		J53	80078578	J, T	441,00											
Altex / Hitachi	100 (HPLC)		R07	80013769		416,00	SD 3651-03		80056123	452,00							
American Hi-Tech	refer to Scientific Systems (SSI)																
Amersham Pharmacia	refer to Pharmacia																
Aminco	DW-2C						SD 3651-03		80056123	452,00							
Aminco	DW-2000						SD 3651-03		80056123	452,00							
Analytik Jena	M500						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	S5 UV/Vis						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	S10 UV/Vis						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 40						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 50						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 50 PLUS						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 200						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 200 PLUS						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 205						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 210						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 210 PLUS						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 250						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord 250 PLUS						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord S 100						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord S 300						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Specord S 600						DO 249 TJ		56070646	J, T	409,00						
Analytik Jena	Spekol 1200	820-60021-0					DO 249 TJ		56070646	J, T	409,00						
Applied Biosystems (ABI)	100S													AB-400	J, L	360,00	
Applied Biosystems (ABI)	120	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00			AB-400	J, L	360,00	
Applied Biosystems (ABI)	120A	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00			AB-400	J, L	360,00	
Applied Biosystems (ABI)	130	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00			AB-400	J, L	360,00	
Applied Biosystems (ABI)	130A	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00			AB-400	J, L	360,00	
Applied Biosystems (ABI)	230A	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00			AB-400	J, L	360,00	
Applied Biosystems (ABI)	270A-HAT (CZE)	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00			AB-400	J, L	360,00	
Applied Biosystems (ABI)	471A	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00						
Applied Biosystems (ABI)	757	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00			AB-400	J, L	360,00	
Applied Biosystems (ABI)	759	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00			AB-400	J, L	360,00	
Applied Biosystems (ABI)	759A	2900-0484					SD 1151-02 TJ		80062563	J, T	460,00			AB-400	J, L	360,00	

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information				Lamp Information									
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Applied Biosystems (ABI)	773	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	775	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00			
Applied Biosystems (ABI)	783	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	783A	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	785A	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	980	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	1000S	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 757	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 759	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 759A	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 773						SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 775						SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 783						SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 783A						SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 785A						SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 980	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00	AB-400	J, L	360,00
Applied Biosystems (ABI)	SF 1000S						SD 1151-02 TJ	80062563	J, T	460,00			
Aqualabo (Secomam)	STAC						DS 253 TJ	45006242	J, T	422,00			
Aqualabo (Secomam)	STAC Chlorine						DS 253 TJ	45006242	J, T	422,00			
Armen Instrument	Spot Prep II 50						SD 1191-01 J	80125095	J	370,00			
Armen Instrument	Spot Prep II 200						SD 1191-01 J	80125095	J	370,00			
Armen Instrument	Spot Prep II 250						SD 1191-01 J	80125095	J	370,00			
Armen Instrument	Spot Skid 3L						SD 1191-01 J	80125095	J	370,00			
Armen Instrument	Spot Skid 5L						SD 1191-01 J	80125095	J	370,00			
Armen Instrument	Spot Skid 6L						SD 1191-01 J	80125095	J	370,00			
ATI	refer to Unicam												
Barspec	all models		R07	80013769		416,00	SD 3651-03	80056123		452,00			
BAS	200		J53	80078578	J, T	441,00							
Beckman	100	585699	R07	80013769	H	416,00	SD 3651-03	80056123		452,00			
Beckman	166	239372	J43	80053989	J, T	626,00							
Beckman	168	538711	J44	80053990	J, T	562,00							
Beckman	DU 600 series	514366	J04	80048278		463,00							
Beckman	DU 620	514366	J04	80048278		463,00							
Beckman	DU 630	514366	J04	80048278		463,00							
Beckman	DU 640	514366	J04	80048278		463,00							
Beckman	DU 650	514366	J04	80048278		463,00							
Beckman	DU 800	514366	J04	80048278		463,00							
Beckman	DU 7000	514366	J04	80048278		463,00							
Beckman	DU 7500	514366	J04	80048278		463,00							
Beckman	P/ACE 5000	144-667					DS 251/05 J	80045466	J	524,00			
Beckman	P/ACE 5010	144-667					DS 251/05 J	80045466	J	524,00			
Beckman	P/ACE 5500	144-667					DS 251/05 J	80045466	J	524,00			
Beckman	P/ACE 5510	144-667					DS 251/05 J	80045466	J	524,00			
Beckman	P/ACE MDQ System	144-667					DS 251/05 J	80045466	J	524,00			
Biocad	60	P5-1085-00-5020	J53	80078578	J, T	441,00							
Biocad	700E	P5-1085-00-5020	J53	80078578	J, T	441,00							
Biocad	Sprint	P5-1085-00-5020	J53	80078578	J, T	441,00							
Biocad	also refer to Perseptive Biosystems												
Biochrom	refer to Pharmacia												
Biorad	206		J53	80078578	J, T	441,00							
Biorad	300		J53	80078578	J, T	441,00							
Biorad	1305		R19	80013778		316,00							
Biorad	1306		R19	80013778		316,00							
Biorad	1503						SD 3651-03	80056123		452,00			
Biorad	1790	930-6106	J53	80078578	J, T	441,00							
Biorad	BioDimensions	930-6106	J53	80078578	J, T	441,00							
Biorad	Biofocus 2000		J53	80078578	J, T	441,00							
Biorad	Biofocus 3000		J53	80078578	J, T	441,00							

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information		Lamp Information											
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Biorad	Biologic Quadtec	7601330						SD 1255-02 TJ	80127318	J, T			378,00
Biorad	CE2000	930-6106	J53	80078578	J, T	441,00							
Biorad	HPE 100	930-6106	J53	80078578	J, T	441,00							
Biorad	MDL 1790	930-6106	J53	80078578	J, T	441,00							
Biorad	MDL 1801	930-6106	J53	80078578	J, T	441,00							
Biorad	Microsampler 100	930-6106	J53	80078578	J, T	441,00							
Bios Analytical	DAD 1100										HP-404	J, L	436,00
Biotage	SP System Flash Chromatography							SD 1191-01 J	80125095	J			370,00
BioTek Instruments	refer to Kontron												
Biotronic	BT3030		R19	80013778		316,00							
Biotronic	BT9510		J53	80078578	J, T	441,00							
Biotronic	BT9520		J53	80078578	J, T	441,00							
Bischoff	DAD3L							SD 1201-01 TJ	80062765	J, T			541,00
Bischoff	Lambda 1010	8250 0200						SD 1201-01 TJ	80062765	J, T			541,00
Büchi	901		R07	80013769		416,00							
Carlo Erba	UV VIS 20		J53	80078578	J, T	441,00							
Carlo Erba	UV VIS 20RS		J53	80078578	J, T	441,00							
Cecil Instruments	1000 series		R18	80078538		397,00							
Cecil Instruments	2000 series		R18	80078538		397,00							
Cecil Instruments	3000 series		R18	80078538		397,00							
Cecil Instruments	5000 series		R18	80078538		397,00							
Cecil Instruments	6000 series		R18	80078538		397,00							
Cecil Instruments	7000 series		R18	80078538		397,00							
Cecil Instruments	8000 series		R18	80078538		397,00							
Cecil Instruments	9000 series		R18	80078538		397,00							
Cecil Instruments	CE 1010	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 1020	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 1070	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 1200	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 1220	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 2021	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 2202	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 2272	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 2292	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 4400	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 5501	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 5502	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 6600	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 6600S	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 6602	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 6603	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 7200	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 7400	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 7460	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 7500	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 7560	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 8020	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 9020	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 9050	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 9080	2202-0142	R18	80078538		397,00							
Cecil Instruments	CE 9090	2202-0142	R18	80078538		397,00							
Cecil Instruments	Series 2	2202-0142	R18	80078538		397,00							
Coleman	55	C055-0505	R07	80013769	H	416,00	SD 3651-03		80056123				452,00
Coleman	124	C055-0505	R07	80013769	H	416,00	SD 3651-03		80056123				452,00
Coleman	550	C055-0505	R07	80013769	H	416,00	SD 3651-03		80056123				452,00
Coleman	552	C055-0505	R07	80013769	H	416,00	SD 3651-03		80056123				452,00
Dionex	2000i	40651	J10	80013703	A	373,00							
Dionex	4000i		J53	80078578	J, T	441,00					LI-401	J, L	323,00

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information						Lamp Information							
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Dionex	AD25	939016T / 6080.2001									PR41065	J, L	452,00
Dionex	CES 1	40651	J10	80013703	A	373,00							
Dionex	CZE	40651	J10	80013703	A	373,00							
Dionex	DSA-1		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Dionex	PDA-100	939016T / 6080.2001									PR41065	J, L	452,00
Dionex	UltiMate (except 3000)						SD 1255-02 TJ	80127318	J, T	378,00			
Dionex	UltiMate DAD-3000	L6999-52 / 6074.1110					XD 6295-01 J	80117195	J, L	452,00			
Dionex	UltiMate DAD-3000RS	L6999-52 / 6074.1110					XD 6295-01 J	80117195	J, L	452,00			
Dionex	UltiMate MWD-3000	L6999-52 / 6074.1110					XD 6295-01 J	80117195	J, L	452,00			
Dionex	UltiMate MWD-3000RS	L6999-52 / 6074.1110					XD 6295-01 J	80117195	J, L	452,00			
Dionex	UltiMate PDA-3000	L9838 / 6080.2001									PR41065	J, L	452,00
Dionex	UltiMate UVD-3000	6073.2070					SD 1255-02 TJ	80127318	J, T	378,00			
Dionex	UltiMate VWD-3000	L6999-52 / 6074.1110					XD 6295-01 J	80117195	J, L	452,00			
Dionex	UltiMate VWD-3100	L6999-52 / 6074.1110					XD 6295-01 J	80117195	J, L	452,00			
Dionex	UltiMate VWD-3400RS	L6999-52 / 6074.1110					XD 6295-01 J	80117195	J, L	452,00			
Dionex	UVD 33	5053.1204									GY-400	J, L	572,00
Dionex	UVD 160	5053.1204					DX 201 RJ	80007875	J, L	467,00	GY-400	J, L	572,00
Dionex	UVD 170S	5053.1204					DX 201 RJ	80007875	J, L	467,00	GY-400	J, L	572,00
Dionex	UVD 170U	5053.1204					DX 201 RJ	80007875	J, L	467,00	GY-400	J, L	572,00
Dionex	UVD 320	5053.1204					DX 201 RJ	80007875	J, L	467,00	GY-400	J, L	572,00
Dionex	UVD 340S	5053.1204					DX 201 RJ	80007875	J, L	467,00	GY-400	J, L	572,00
Dionex	UVD 340U	5053.1204					DX 201 RJ	80007875	J, L	467,00	GY-400	J, L	572,00
Dionex	VDM-1		J53	80078578	J, T	441,00							
Dupont	825						SD 3651-03	80056123		452,00			
Dupont	850						SD 3651-03	80056123		452,00			
Dynamax	UV-C		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Dynamax	UV-D11		J53	80078578	J, T	441,00							
Dynamax	UV-M		J53	80078578	J, T	441,00					LI-401	J, L	323,00
ECOM	Flash 06	22590001					SD 1191-01 J	80125095	J	370,00			
ECOM	Flash DAD Series	22590005					SD 1191-01 J	80125095	J	370,00			
ECOM	LCD 2073A	22590001					SD 1191-01 J	80125095	J	370,00			
ECOM	LCD 2083.2	22590000					SD 1291-01 J	80087724	J	382,00			
ECOM	LCD 2083.2 CE	22590000					SD 1291-01 J	80087724	J	382,00			
ECOM	LCD 2083.61	22590000					SD 1291-01 J	80087724	J	382,00			
ECOM	LCD 2084.2	22590000					SD 1291-01 J	80087724	J	382,00			
ECOM	LCD 2084.2 CE	22590000					SD 1291-01 J	80087724	J	382,00			
ECOM	Opal	22590001					SD 1191-01 J	80125095	J	370,00			
ECOM	Sapphire	22590001					SD 1191-01 J	80125095	J	370,00			
ECOM	Topaz	22590001					SD 1191-01 J	80125095	J	370,00			
EM Science	all models						SD 3651-04 J	80056072	J	395,00			
Erma	ERC 7210						SD 3651-03	80056123		452,00			
Erma	ERC 7211						SD 3651-03	80056123		452,00			
Filmtrics	F20-UV (1998 and earlier)	LAMP-D2-DH2000					XD 5665-11 J	80087387	J, L	678,00			
Filmtrics	LS-DT2	LAMP-D2-LS-DT2					XD 5685-01 J	80058872	J, L	636,00			
GAT	LCD 500		J53	80078578	J, T	441,00							
GAT	LCD 502		J53	80078578	J, T	441,00							
GAT	LCD 503		J53	80078578	J, T	441,00							
GAT	LCD 550		J53	80078578	J, T	441,00							
GAT	PHD 601		J53	80078578	J, T	441,00							
GBC	LC1205		J53	80078578	J, T	441,00					LI-401	J, L	323,00
GBC	LC1206		J53	80078578	J, T	441,00					LI-401	J, L	323,00
GBC	LC1210		J53	80078578	J, T	441,00					LI-401	J, L	323,00
GE Healthcare	Akta Prime										GE-400	L	589,00
Gilford	Response		R07	80013769		416,00	SD 3651-03	80056123		452,00			
Gilson	115	100326					PR38080	80024223	J	497,00			
Gilson	116	100326					PR38080	80024223	J	497,00			
Gilson	117	100326					PR38080	80024223	J	497,00			
Gilson	118	100326					PR38080	80024223	J	497,00			

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information				Lamp Information									
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Gilson	119	100326					PR38080	80024223	J	497,00			
Gilson	151	100326					PR38080	80024223	J	497,00			
Gilson	152	100326					PR38080	80024223	J	497,00			
Gilson	153	100326					PR38080	80024223	J	497,00			
Gilson	155	100326					PR38080	80024223	J	497,00			
Gilson	156	100326					PR38080	80024223	J	497,00			
Gilson	160 DAD	1G041402					PR38080	80024223	J	497,00			
Gilson	170 DAD	18011003					PR38080	80024223	J	497,00			
Gynkotek	refer to Dionex												
Hach Lange	DR 5000	A23792 / 2986000					XD 6251-01 J	80093632	J, L	414,00			
Hach Lange	DR 6000	A23792 / 2986000					XD 6251-01 J	80093632	J, L	414,00			
Hewlett Packard	refer to Agilent Technologies												
Hitachi	100-10	982-1035 or 799-9991	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Hitachi	100-30	982-1035 or 799-9991	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Hitachi	100-40	982-1035 or 799-9991	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Hitachi	100-50	982-1035 or 799-9991	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Hitachi	100-60	982-1035 or 799-9991	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Hitachi	124	982-1035 or 799-9991	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Hitachi	150-20	J 38-4017					SD 3651-04 J	80056072	J	395,00			
Hitachi	181	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	200	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	220	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	220A	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	300	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	330	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	330LC	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	340	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	655A22	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	5410 Chromaster	892-2550					SD 3651-04 J	80056072	J	395,00			
Hitachi	5420 Chromaster	892-2550					SD 3651-04 J	80056072	J	395,00			
Hitachi	5430 Chromaster	892-2550					SD 3651-04 J	80056072	J	395,00			
Hitachi	L-2000	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	L-2400 (LCE)	890-2430									0342-098	J, L	406,00
Hitachi	L-2420 (LCE)	890-2430									0342-098	J, L	406,00
Hitachi	L-2450	890-2430									0342-098	J, L	406,00
Hitachi	L-2455 DAD	890-2430									0342-098	J, L	406,00
Hitachi	L-2500	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-3000	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-3200	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	L-4000	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-4200	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-4225	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-4250	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-4500	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-4520	371140005									0342-098	J, L	406,00
Hitachi	L-4720	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-6000	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-7000	890-2430									MH-404	J, L	323,00
Hitachi	L-7400	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-7420	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-7450	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-7450A	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	L-7455	885-3570 or 890-2430									0342-098	J, L	406,00
Hitachi	LaChrom										0342-098	J, L	406,00
Hitachi	U-1000	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-1100	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-1500	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-1800	122-2300					SD 3651-07 J	80063164	J	395,00			

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information				Lamp Information									
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Hitachi	U-2000	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-2001	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-2800	122-2300					SD 3651-07 J	80063164	J	395,00			
Hitachi	U-2810	122-2300					SD 3651-07 J	80063164	J	395,00			
Hitachi	U-2900	122-2300					SD 3651-07 J	80063164	J	395,00			
Hitachi	U-3010	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-3200	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-3210	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-3310	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-3410	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-3501	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-4000	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-4001	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-4010	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	U-4100	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi	UL-1100	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi / Perkin Elmer	124	239-0354					SD 3651-04 J	80056072	J	395,00			
Hitachi / Perkin Elmer	139		R07	80013769		416,00	SD 3651-03	80056123		452,00			
Hitachi / Perkin Elmer	200	239-0354					SD 3651-04 J	80056072	J	395,00			
Hitachi / Perkin Elmer	320	239-0354					SD 3651-04 J	80056072	J	395,00			
Hitachi / Perkin Elmer	329	239-0354					SD 3651-04 J	80056072	J	395,00			
Hitachi / Perkin Elmer	330	885-1000					SD 3651-04 J	80056072	J	395,00			
Hitachi / Perkin Elmer	340	239-0354					SD 3651-04 J	80056072	J	395,00			
Hitachi / Perkin Elmer	556	239-0354					SD 3651-04 J	80056072	J	395,00			
Hitachi / Perkin Elmer	557	239-0354					SD 3651-04 J	80056072	J	395,00			
IBM	Spectrophotometers		R19	80013778		316,00							
ICI	LC1200		R07	80013769	H	416,00	SD 3651-03	80056123		452,00			
ICI	LC1205		J53	80078578	J, T	441,00							
ICI	LC1210		J53	80078578	J, T	441,00							
IFCO	S500		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Interchim	Puriflash 215 (GO9630)						SD 1191-01 J	80125095	J	370,00			
Interchim	Puriflash 430 (GO9630)						SD 1191-01 J	80125095	J	370,00			
Interchim	Puriflash 450 (GC9630)						SD 1191-01 J	80125095	J	370,00			
Interchim	Puriflash 4100 (GO9630)						SD 1191-01 J	80125095	J	370,00			
Interchim	Puriflash 4250-40 (GO9630)						SD 1191-01 J	80125095	J	370,00			
Interchim	Puriflash 4250-250 (GO9630)						SD 1191-01 J	80125095	J	370,00			
Jasco	975 (B series)	L6300-51									JV-401	J, L	323,00
Jasco	975 (C series)	L6300-51									JV-401	J, L	323,00
Jasco	1575	5330-0091									JV-401	J, L	323,00
Jasco	2050		R19	80013778		316,00							
Jasco	2075	5330-0091									JV-401	J, L	323,00
Jasco	Multi 320		R19	80013778		316,00							
Jasco	Uvidec 4D		R07	80013769		416,00	SD 3651-03	80056123		452,00			
Jasco	Uvidec 100 III		R19	80013778		316,00							
Jasco	Uvidec 100 IV		R19	80013778		316,00							
Jasco	Uvidec 100 V		R19	80013778		316,00							
Jasco	Uvidec 100 VI		R19	80013778		316,00							
Jasco	Uvidec 430		R19	80013778		316,00							
Jasco	Uvidec 610		R19	80013778		316,00							
Jasco	Uvidec 650		R19	80013778		316,00							
Jasco	V-530										WL24443 A	J	274,00
Jasco	V-550										WL24443 A	J	274,00
Jasco	V-630										WL24443 A	J	274,00
Jasco	V-630Bio										WL24443 A	J	274,00
Jasco	V-650										WL24443 A	J	274,00
Jasco	V-660										WL24443 A	J	274,00
Jasco	V-670										WL24443 A	J	274,00
Kittiwake	Procal 5000 Ultra-Violet Analyser		Y06	80017708		717,00							

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information						Lamp Information							
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Knauer	DAD K-2800	A4447					DX 241/05 TJ	80004738	J, L, T	488,00			
Knauer	Smartline PDA Detector 2800	A4447					DX 241/05 TJ	80004738	J, L, T	488,00			
Knauer	Smartline PDA Detector 2850	A4447					DX 241/05 TJ	80004738	J, L, T	488,00			
Knauer	Smartline UV Detector 2500	A4071-1					SD 1255-02 TJ	80127318	J, T	378,00			
Knauer	Smartline UV Detector 2600	A4071-1					SD 1255-02 TJ	80127318	J, T	378,00			
Knauer	WellChrom K-2000	A4071-1					SD 1255-02 TJ	80127318	J, T	378,00	KA-401	J, L	323,00
Knauer	WellChrom K-2001	A4071-1					SD 1255-02 TJ	80127318	J, T	378,00	KA-401	J, L	323,00
Knauer	WellChrom K-2500	A4071-1					SD 1255-02 TJ	80127318	J, T	378,00	KA-401	J, L	323,00
Knauer	WellChrom K-2501	A4071-1					SD 1255-02 TJ	80127318	J, T	378,00	KA-401	J, L	323,00
Knauer	WellChrom K-2600	A4071-1					SD 1255-02 TJ	80127318	J, T	378,00	KA-401	J, L	323,00
Knauer	WellChrom K-2700	A4071-1					SD 1255-02 TJ	80127318	J, T	378,00			
Kontron	332	91-00541A					DS 222/05 TJ	80006939	J, T	447,00			
Kontron	335	91-91494					DX 222/05 TJ	80005899	J, L, T	488,00			
Kontron	430	91-00541A					DS 222/05 TJ	80006939	J, T	447,00			
Kontron	430A	91-91494					DX 222/05 TJ	80005899	J, L, T	488,00			
Kontron	432	91-00541A					DS 222/05 TJ	80006939	J, T	447,00			
Kontron	433	91-00541A					DS 222/05 TJ	80006939	J, T	447,00			
Kontron	535	91-91494					DX 222/05 TJ	80005899	J, L, T	488,00			
Kontron	535 DAD	91-91494					DX 222/05 TJ	80005899	J, L, T	488,00			
Kontron	535 VVD	91-91494					DX 222/05 TJ	80005899	J, L, T	488,00			
Kontron	540+	54-02007					DX 241/05 TJ	80004738	J, L, T	488,00			
Kontron	540 DAD	54-02007					DX 241/05 TJ	80004738	J, L, T	488,00			
Kontron	545V	54-02007					DX 241/05 TJ	80004738	J, L, T	488,00			
Kontron	720 LC						DO 905 T	45002279	T, mBS	373,00			
Kontron	725						SD 3651-03	80056123		452,00			
Kontron	725 LC	KO 3017	R07	80013769		416,00							
Kontron	730 LC	KO 0302004	R19	80013778		316,00							
Kontron	735		R19	80013778		316,00							
Kontron	735 LC	9600 2746	R19	80013778		316,00							
Kontron	Uvikon 430						DS 222/05 TJ	80006939	J, T	447,00			
Kontron	Uvikon 432						DS 222/05 TJ	80006939	J, T	447,00			
Kontron	Uvikon 720						DO 905 T	45002279	T, mBS	373,00			
Kontron	Uvikon 720 LC						DO 905 T	45002279	T, mBS	373,00			
Kontron	Uvikon 722		R19	80013778		316,00							
Kontron	Uvikon 725		R07	80013769		416,00	SD 3651-03	80056123		452,00			
Kontron	Uvikon 730		R19	80013778		316,00							
Kontron	Uvikon 900	90-00782ST					SD 3611-02 TJ	80052431	J, T	496,00			
Kontron	Uvikon 922	90-00782ST					SD 3611-02 TJ	80052431	J, T	496,00			
Kontron	Uvikon 923	90-00782ST					SD 3611-02 TJ	80052431	J, T	496,00			
Kontron	Uvikon 930	90-00782ST					SD 3611-02 TJ	80052431	J, T	496,00			
Kontron	Uvikon 931	90-00782ST					SD 3611-02 TJ	80052431	J, T	496,00			
Kontron	Uvikon 932	90-00782ST					SD 3611-02 TJ	80052431	J, T	496,00			
Kontron	Uvikon 933	90-00782ST					SD 3611-02 TJ	80052431	J, T	496,00			
Kontron	Uvikon 940	90-00782ST					SD 3611-02 TJ	80052431	J, T	496,00			
Kratos	refer to Applied Biosystems												
Kronlab (YMC)	KLA.2001	KLA.4071					SD 1255-02 TJ	80127318	J, T	378,00			
Kronlab (YMC)	KLA.2501	KLA.4071					SD 1255-02 TJ	80127318	J, T	378,00			
Labomatic	Labocord-3000						XD 3465-01 TJ	80127418	J, L, T	444,00			
Labomatic	Labocord-5000						XD 3465-01 TJ	80127418	J, L, T	444,00			
LC Packings	refer to Dionex												
Linear	200		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Linear	201		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Linear	203		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Linear	204		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Linear	205		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Linear	206		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Linear	206PH		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Linear	Focus		J53	80078578	J, T	441,00					LI-401	J, L	323,00

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information							Lamp Information						
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Linear	Focus 2000										LI-401	J, L	323,00
Linear	Focus 2000 HR	9551-0202	J53	80078578	J, T	441,00					LI-401	J, L	323,00
Linear	UV 100										LI-401	J, L	323,00
Linear	UV 200		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Linear	UV 203										LI-401	J, L	323,00
Linear	UV 204										LI-401	J, L	323,00
LKB Pharmacia	4050		R07	80013769	E	416,00	SD 3651-03	80056123		452,00			
LKB Pharmacia	4054		R07	80013769	E	416,00	SD 3651-03	80056123		452,00			
LKB Pharmacia	TDS (Ser# > 205)		R07	80013769	E	416,00	SD 3651-03	80056123		452,00			
Merck / Hitachi	refer to Hitachi												
Metertech	SP-8001						SD 1155-01 J	80082959		457,00			
Micronal	B582										WL24443 A	J	274,00
Molecular Devices	FilterMax F-series	S PD 400 001									PR38033	J	753,00
Nicolet	8745		L06	80013738		320,00							
Nicolet	Evolution 100						SD 1251-03 J	800135828	J	374,00			
Nicolet	Evolution 500						SD 1251-03 J	800135828	J	374,00			
Nicolet	Evolution 600						SD 1251-03 J	800135828	J	374,00			
Perkin Elmer	55	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	550	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	551	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	552	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	552A	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	553	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	559	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	559A	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	575	C0550505					SD 3651-03	80056123		452,00			
Perkin Elmer	785A	2900-0484					SD 1151-02 TJ	80062563	J, T	460,00			
Perkin Elmer	Integral 2000	N2351285	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	Integral 4000	N2351285	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	Lambda 1	C0550505					SD 3651-03	80056123		452,00			
Perkin Elmer	Lambda 2	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 2S	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 3	C0550505					SD 3651-03	80056123		452,00			
Perkin Elmer	Lambda 3A	C0550505					SD 3651-03	80056123		452,00			
Perkin Elmer	Lambda 3B	C0550505					SD 3651-03	80056123		452,00			
Perkin Elmer	Lambda 4	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 5	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 7	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 9	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 10	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 11	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 12	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 14	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 14 P	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 15	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 16	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 17	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 18	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 19	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 20	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 25	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 35	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 40	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 40 P	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 45	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 650	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	455,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 800	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	455,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 850	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information		Lamp Information											
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Perkin Elmer	Lambda 900	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 950	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda Bio 19	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda Bio 20	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda Bio 40	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda 365	N4104016					XD 3495-02 J	80138209	J	480,00			
Perkin Elmer	Lambda 465	N4104016					XD 3495-02 J	80138209	J	480,00			
Perkin Elmer	Lambda Bio 800	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda Bio 900	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Lambda EZ150	C0550505 or C6957529					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer	Lambda EZ201	C0550505 or C6957529					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer	Lambda EZ210	C0550505 or C6957529					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer	LC-55	02711340	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-55B	02711340	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-65B	02711340	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-65T	02711340	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-75	02711340	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-85	02711340 or N2351285	R07	80013769		416,00							
Perkin Elmer	LC-85B	02711340	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-90	0271-2224											
Perkin Elmer	LC-95	02711340	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-135C	N2351285	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-141MC		R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-235C	N2351285	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-241MC	02712224	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer	LC-295	02712266					SD 1151-02 TJ	80062563	J, T	460,00			
Perkin Elmer	LC-480 Photo Diode Array	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	LC Series 200 Diode Array Detector, Phase 2	N2922046									PR38058	J	477,00
Perkin Elmer	M55	B0160917 or C055-0505					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	M550	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	M551	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	M552 (Ser# > 2030)	102014					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	M552A	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	M553	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	M559	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	M559A	B0160917 or L6022728					SD 3651-06 TJ	80055323	J, T	480,00	replaces DO 650 TJ (80007103)		
Perkin Elmer	Series 200 DAD, Phase I	N2922010					DS 272 TJ	45002219	J	490,00			
Perkin Elmer	Series 200 DAD, Phase II	N2922046					DS 272 TJ	45002219	J	490,00			
Perkin Elmer	Series 200 UV/Vis	N2900484					SD 1151-02 TJ	80062563	J, T	460,00			
Perkin Elmer	Series LC 200 DAD	N2922010					DS 272 TJ	45002219	J, T	490,00			
Perkin Elmer / ABI	refer to Applied Biosystems												
Perkin Elmer / Hitachi	124	239-0354					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer / Hitachi	139	98201035	R07	80013769		416,00	SD 3651-03	80056123		452,00			
Perkin Elmer / Hitachi	200	239-0354					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer / Hitachi	320	239-0354					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer / Hitachi	329	239-0354					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer / Hitachi	330	239-0354					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer / Hitachi	340	239-0354					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer / Hitachi	556	239-0354					SD 3651-04 J	80056072	J	395,00			
Perkin Elmer / Hitachi	557	239-0354					SD 3651-04 J	80056072	J	395,00			
Perseptive Biosystems	Biocad 20		J53	80078578	J, T	441,00							
Perseptive Biosystems	Biocad 60		J53	80078578	J, T	441,00							
Perseptive Biosystems	BioCAD/SPRINT system for Perfusion Chromatography		J53	80078578	J, T	441,00							
Pharmacia	refer to LKB Pharmacia												
Phillips	refer to Unicam												
Pye Unicam	refer to Unicam												
Scientific Systems Inc (SSI)	200 UV		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Scientific Systems Inc (SSI)	201		J53	80078578	J, T	441,00					LI-401	J, L	323,00

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information				Lamp Information										
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price	
Scientific Systems Inc (SSI)	500		J53	80078578	J, T	441,00					LI-401	J, L	323,00	
Scientific Systems Inc (SSI)	506 UV/Vis		J53	80078578	J, T	441,00					LI-401	J, L	323,00	
Scientific Systems Inc (SSI)	525		J53	80078578	J, T	441,00					LI-401	J, L	323,00	
Scientific Systems Inc (SSI)	UV 200		J53	80078578	J, T	441,00					LI-401	J, L	323,00	
Scinco	S-3100							XD 3495-03 J	80139864	J, L			392,00	
Scinco	S-4100							XD 3495-03 J	80139864	J, L			392,00	
Severn Analytical	6500										LI-401	J, L	323,00	
Severn Analytical	6503										LI-401	J, L	323,00	
Severn Analytical	6504										LI-401	J, L	323,00	
Severn Analytical	6508										LI-401	J, L	323,00	
Severn Analytical	SA6500		J53	80078578	J, T	441,00								
Severn Analytical	SA6503		J53	80078578	J, T	441,00								
Severn Analytical	SA6504		J53	80078578	J, T	441,00								
Severn Analytical	SA6508 (Rapidsan)		J53	80078578	J, T	441,00								
Severn Analytical	SA6510							SD 1151-02 TJ	80062563	J, T			460,00	
Shimadzu	AA-670	062-65055-05									WL24443 A	J	274,00	
Shimadzu	AA-670G	062-65055-05									WL24443 A	J	274,00	
Shimadzu	AA-6200	062-65055-05									WL24443 A	J	274,00	
Shimadzu	AA-6300	062-65055-05									WL24443 A	J	274,00	
Shimadzu	AA-6500	062-65055-05									WL24443 A	J	274,00	
Shimadzu	AA-6701	062-65055-05									WL24443 A	J	274,00	
Shimadzu	AA-6800	062-65055-05									WL24443 A	J	274,00	
Shimadzu	CS-920	062-65055-05									WL24443 A	J	274,00	
Shimadzu	CS-930	062-65055-05									WL24443 A	J	274,00	
Shimadzu	CS-9000	062-65055-05									WL24443 A	J	274,00	
Shimadzu	D300L	062-65055-05									WL24443 A	J	274,00	
Shimadzu	LC-3A										SH-400	L	323,00	
Shimadzu	LC-4A										SH-400	L	323,00	
Shimadzu	LC-6A										SH-400	L	323,00	
Shimadzu	LC-2010	228-37401-00						XD 6295-07 J	80137686	J, L			367,00	
Shimadzu	LC-2010A	228-37401-00						XD 6295-07 J	80137686	J, L			367,00	
Shimadzu	LC-2010C	228-37401-00						XD 6295-07 J	80137686	J, L			367,00	
Shimadzu	MPS-2000	062-65055-05									WL24443 A	J	274,00	
Shimadzu	SP-4										SH-400	L	323,00	
Shimadzu	SPD-2A	670-14501-02						XD 3441-01	80139739	L	425,00	SH-400	L	323,00
Shimadzu	SPD-2AV										SH-400	L	323,00	
Shimadzu	SPD-6A	062-65056-03						XD 3441-01	80139739	L	425,00	SH-400	L	323,00
Shimadzu	SPD-6AV	062-65056-03						XD 3441-01	80139739	L	425,00	SH-400	L	323,00
Shimadzu	SPD-7A	062-65056-03						XD 3441-01	80139739	L	425,00	SH-400	L	323,00
Shimadzu	SPD-7AV	062-65056-03						XD 3441-01	80139739	L	425,00	SH-400	L	323,00
Shimadzu	SPD-10	228-34016-02						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	SPD-10A	228-34016-02						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	SPD-10AS	228-34016-02						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	SPD-10AV	228-34016-02						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	SPD-10AVP	228-34016-02						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	SPD-10AVVP	228-34016-02						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	SPD-20A	228-34016-02						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	SPD-20AV	228-34016-02						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	SPD-M10A	670-14604-00						DX 243/05 J	80004249	J, L			555,00	
Shimadzu	SPD-M10AV	670-14604-00						DX 243/05 J	80004249	J, L			555,00	
Shimadzu	SPD-M10AVP	228-34016-00						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	SPD-M20A	228-34016-00						DX 250/05 J	80003115	J, L			392,00	
Shimadzu	UV-120	200-75503-01									WL24443 A	J	274,00	
Shimadzu	UV-120-01	200-75503-01									WL24443 A	J	274,00	
Shimadzu	UV-120-02	200-75503-01									WL24443 A	J	274,00	
Shimadzu	UV-140	200-75503-01									WL24443 A	J	274,00	
Shimadzu	UV-150	200-75503-01									WL24443 A	J	274,00	
Shimadzu	UV-160	200-75503-01									WL24443 A	J	274,00	

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information				Lamp Information									
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Shimadzu	UV-160A	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-180	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-190	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-200	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-202	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-210	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-210A	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-240	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-250	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-260	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-265	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-350	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-360	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-365	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-370	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-625	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-1200	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-1201	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-1600	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-1700	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-1800	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-2100	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-2101	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-2400	062-65055-05									WL24443 A	J	274,00
Shimadzu	UV-2401 PC	062-65055-05									WL24443 A	J	274,00
Shimadzu	UV-2450	062-65055-05									WL24443 A	J	274,00
Shimadzu	UV-2500	062-65055-05									WL24443 A	J	274,00
Shimadzu	UV-2501 PC	062-65055-05									WL24443 A	J	274,00
Shimadzu	UV-2550	062-65055-05									WL24443 A	J	274,00
Shimadzu	UV-2600	200-75503-02									WL24443 A	J	274,00
Shimadzu	UV-2700	200-75503-02									WL24443 A	J	274,00
Shimadzu	UV-3000	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-3001	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-3100	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-3101	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-3600	200-75503-01									WL24443 A	J	274,00
Shimadzu	UV-3700	200-75503-01									WL24443 A	J	274,00
Shimadzu	Uvmini-1240	200-75503-01									WL24443 A	J	274,00
Speck Analytical	SA 203		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Speck Analytical	SA 302		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	100		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	200		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	CE 100		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	CE 500		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	CE 1000		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	CE 2000		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	Focus Detector	9551-0202	J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	SP 100		J53	80078578	J, T	441,00							
Spectra Physics	SP 200		J53	80078578	J, T	441,00							
Spectra Physics	SP 770							SD 1151-02 TJ	80062563	J, T	460,00		
Spectra Physics	SP 970							SD 1151-02 TJ	80062563	J, T	460,00		
Spectra Physics	SP 8480							SD 1151-02 TJ	80062563	J, T	460,00		
Spectra Physics	Spectrachrom 100		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	Spectrachrom 200		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	UV100		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	UV150		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	UV200		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	UV1000		J53	80078578	J, T	441,00					LI-401	J, L	323,00

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information				Lamp Information									
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Spectra Physics	UV2000		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Spectra Physics	UV3000										LI-401	J, L	323,00
Spectra Physics	UV6000LP	108052						SD 1255-01 TJ	80122679	J, T	512,00		
Spectra Physics	also refer to Thermo												
Spectro BioNova	Celeno DAD							DX 241/05 TJ	80004738	J, L, T	488,00		
Spectronic Instruments	refer to Thermo												
Suntek Science Instruments	refer to Scientific Systems Inc												
Therma-Wave	Series 3	29-013586						D6255-01J	80133979	J	424,00		
Therma-Wave	Series 5	29-013586						D6255-01J	80133979	J	424,00		
Thermo	Aquamate							SD 1251-03 J	800135828	J	374,00		
Thermo	Biomate							SD 1251-03 J	800135828	J	374,00		
Thermo	Evolution 600							SD 1251-03 J	800135828	J	374,00		
Thermo	Genesis 8	9423-UV9-002E / 9004B						SD 1251-03 J	800135828	J	374,00		
Thermo	iCE 3300 AAS	9423 420 30004						SD 1241-01 J	80119072	J	320,00		
Thermo	M Series (older series)	9423 420 20004									PR38002	J	428,00
Thermo	M Serie (newer series)	9423 420 30004						SD 1241-01 J	80119072	J	320,00		
Thermo	S Series	9423 420 30004						SD 1241-01 J	80119072	J	320,00		
Thermo	SC 100		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	SC 200		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	Solaar M5	9423 420 20004									PR38002	J	428,00
Thermo	Solaar M6	9423 420 20004									PR38002	J	428,00
Thermo	Solaar S4	9423 420 30004						SD 1241-01 J	80119072	J	320,00		
Thermo	Spectra 100	9551-0023	J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	Spectrachrom 100		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	Spectrachrom 200		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	Surveyor Plus	108052						SD 1255-01 TJ	80122679	J, T	512,00		
Thermo	UV1							SD 1251-03 J	800135828	J	374,00		
Thermo	UV1-100							SD 1251-03 J	800135828	J	374,00		
Thermo	UV2							SD 1251-03 J	800135828	J	374,00		
Thermo	UV3							SD 1251-03 J	800135828	J	374,00		
Thermo	UV4							SD 1251-03 J	800135828	J	374,00		
Thermo	UV100	9551-0023	J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	UV120										LI-401	J, L	323,00
Thermo	UV150	9551-0023	J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	UV200	9551-0023	J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	UV1000	9551-0023	J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	UV2000	9551-0023	J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	UV3000	9551-0023	J53	80078578	J, T	441,00					LI-401	J, L	323,00
Thermo	UV3000HR	9551-0202 - A4903-010	J53	80078578	J, T	441,00							
Thermo	UV6000LP	108052						SD 1255-01 TJ	80122679	J, T	512,00		
Thermo	also refer to Spectra Physics and Unicam												
Thermo Elemental	refer to Thermo												
Thermo Finnigan	refer to Thermo												
Thermo Quest	refer to Thermo												
Thermo Separation	refer to Thermo												
Thermo Spectronic	refer to Thermo												
TJA Solutions	refer to Thermo and Unicam												
Unicam	929 Atomic Absorption Spectrometer	9423 393 46251									PR18169		428,00
Unicam	939 Atomic Absorption Spectrometer	9423 393 46251									PR18169		428,00
Unicam	959 Atomic Absorption Spectrometer	9423 393 46251									PR18169		428,00
Unicam	969 Atomic Absorption Spectrometer	9423 393 46251									PR18169		428,00
Unicam	4225		J53	80078578	J, T	441,00					LI-401	J, L	323,00
Unicam	Evolution 600							SD 1251-03 J	800135828	J	374,00		
Unicam	Helios II	9423-UV9-002E / 9004B						SD 1251-03 J	800135828	J	374,00		
Unicam	Helios Alpha	9423-UV9-002E / 9004B						SD 1251-03 J	800135828	J	374,00		
Unicam	Helios Beta	9423-UV9-002E / 9004B						SD 1251-03 J	800135828	J	374,00		
Unicam	Helios Delta	9423-UV9-002E / 9004B						SD 1251-03 J	800135828	J	374,00		
Unicam	Helios Epsilon	9423-UV9-002E / 9004B						SD 1251-03 J	800135828	J	374,00		

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information							Lamp Information						
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
Unicam	Helios Gamma	9423-UV9-002E / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	Helios Omega	9423-UV9-002E / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	Helios Zeta	9423-UV9-002E / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	SP8450		L06	80013738	B	320,00							
Unicam	SP8700		L06	80013738	B	320,00							
Unicam	SP8710		L06	80013738	B	320,00							
Unicam	SP8715		L06	80013738	B	320,00							
Unicam	SP8720		L06	80013738	B	320,00							
Unicam	SP8730		L06	80013738	B	320,00							
Unicam	SP8740	4013-163-75402	L06	80013738	B	320,00							
Unicam	SP8745		L06	80013738	B	320,00							
Unicam	SP8750		L06	80013738	B	320,00							
Unicam	UV1	9423U / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	UV2	9423U / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	UV3	9423U / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	UV4	9423U / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	UV300	9423U / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	UV500	9423U / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	UV550	9423U / 9004B					SD	1251-03 J	800135828	J			374,00
Unicam	also refer to Thermo												
Varian	2050	03-915303-99	R19	80013778		316,00							
Varian	9050	391615691					XD	3445-02 TJ	80003865	J, L, T			571,00
Varian	ProStar 310	391615691					XD	3445-02 TJ	80003865	J, L, T			571,00
Varian	ProStar 325	1107 15400					DX	230/05 J	80051354	J, L			541,00
Varian	ProStar 330 DAD	393570502									0342-098	J, L	406,00
Varian	ProStar 335	1107 15400					DX	230/05 J	80051354	J, L			541,00
Varian	ProStar 340	R000088515	J53	80078578	J, T	441,00							
Varian	ProStar 345	R000088515	J53	80078578	J, T	441,00							
Varian	UV-50	391615691					XD	3445-02 TJ	80003865	J, L, T			571,00
Varian	UV-100	391615691					XD	3445-02 TJ	80003865	J, L, T			571,00
Varian	UV-200	391615691					XD	3445-02 TJ	80003865	J, L, T			571,00
Varian	Vista 5000	391615691					XD	3445-02 TJ	80003865	J, L, T			571,00
Varian	Vista 5500	391615691					XD	3445-02 TJ	80003865	J, L, T			571,00
Waters	480LC	WAT099499					DO	903 T	56069329	T			399,00
Waters	481	WAT099499					DO	903 T	56069329	T			399,00
Waters	481 Lambda Max	WAT099499					DO	903 T	56069329	T			399,00
Waters	481LC	WAT099499					DO	903 T	56069329	T			399,00
Waters	486	WAT080678									PR41056	J, L, T	443,00
Waters	486 LC Module 1	WAT080678									WA-402	J, L	426,00
Waters	487	WAS081142					XD	6295-03 J	80119133	J, L			455,00
Waters	996 DAD	WAT052586 (f. Longlife)					DX	247/05 TJ	80006575	J, L, T			488,00
Waters	2487 Dual Wavelength	WAS081142					XD	6295-03 J	80119133	J, L			455,00
Waters	2488 Dual Wavelength	WAS081142					XD	6295-03 J	80119133	J, L			455,00
Waters	2996 PDA	WAT052586					DX	247/05 TJ	80006575	J, L, T			488,00
Waters	ACQUITY 2996 PDA Detector	WAT052586					DX	247/05 TJ	80006575	J, L, T			488,00
Waters	Alliance 2487	WAS081142					XD	6295-03 J	80119133	J, L			455,00
Waters	Alliance 2488	WAS081142					XD	6295-03 J	80119133	J, L			455,00
Waters	UV2487	WAS081142					XD	6295-03 J	80119133	J, L			455,00
Zeiss (Jena)	refer to Analytik Jena												

Deuterium Lamp Application List

All prices are in Euros per piece and are subject to change without further notice.
All brand names and registered trademarks are the property of their respective owners.

Instrument Information		Lamp Information											
Instrument Manufacturer	Model	OEM	Type	Part No.	Code	Price	Type	Part No.	Code	Price	Type/Part No.	Code	Price
A - use clip from old lamp B - use bracket from old lamp C - requires prealigning D - use connector from old lamp E - resolder circuit board from old lamp onto new lamp F - check previous lamp G - bare lamp H - check terminals J - prealigned L - longlife (2000 hrs) N - exchange only T - with timer X - pos. adapter req. mBS - with Burndy connector mGK - with spade lugs mK - with cooling jacket oK - without cooling jacket mMS - with Molex connector nR - without RFID													

Heraeus Fiberlight®

All prices are in Euros per piece and are subject to change without further notice.
 All brand names and registered trademarks are the property of their respective owners.

Description	Part No.	Price
FiberLight Modules		
FiberLight D2 Basic DTM 6/10 - aperture size: 1.0 mm, quartz glass, focused beam, for 400 - 600 µm optical fiber, 200 - 1100 nm	80007444	775,00
FiberLight D2 Basic DTM 6/50 - aperture size: 0.5 mm, quartz glass, focused beam, for 200 - 600 µm optical fiber, 200 - 1100 nm	45006281	775,00
FiberLight D2 Basic DTM 6/10S - aperture size: 1.0 mm, synth. quartz glass, focused beam, for 400 - 600 µm optical fiber, 185 - 1100 nm	80000755	775,00
FiberLight D2 Basic DTM 6/50S - aperture size: 0.5 mm, synth. quartz glass, focused beam, for 200 - 600 µm optical fiber, 185 - 1100 nm	80001017	775,00
FiberLight D2 Basic DTM 6/11 - aperture size: 1.0 mm, quartz glass, collimated beam (quasi parallel light)	80003030	775,00
FiberLight D2 Basic DTM 6/11S - aperture size: 1.0 mm, synth. quartz glass, collimated beam (quasi parallel light)	80048581	775,00
FiberLight D2 Basic DTM 6/15 - aperture size: 1.0 mm, quartz glass, focused beam, for 400 - 600 µm optical fiber, 200 - 1100 nm	80010770	775,00
FiberLight D2 Compact DTM 6/17 - aperture size: 1.0 mm, quartz glass, focused beam, for 400 - 600 µm optical fiber, 200 - 400 nm	80012065	775,00
FiberLight D2 HighPower DTM 10/50S - aperture size: 0.5 mm, synth. quartz glass, focused beam, for 200 - 600 µm optical fiber, 185 - 1100 nm	80043681	1.133,00
Replacement Lamp Units		
FiberLight D2 Replacement Lamp DTL 6/10, for DTM 6/10, DTM 6/11	45006253	419,00
FiberLight D2 Replacement Lamp DTL 6/50, for DTM 6/50	45006266	419,00
FiberLight D2 Replacement Lamp DTL 6/10S, for DTM 6/10S, DTM 6/11S	80000756	419,00
FiberLight D2 Replacement Lamp DTL 6/50S, for DTM 6/50S	80001018	419,00
FiberLight D2 Replacement Lamp DTL 6/15, for DTM 6/15	80010771	419,00
FiberLight D2 Replacement Lamp DTL 10/50S, for DTM 10/50S	80088565	568,00