

## pH Indicator Test Papers

### Advantages of Test Papers

- Determination of pH value easily and inexpensively, without instrumentation
- Color changes instantaneously, compare to chart for determination
- Two convenient formats:
  - Roll, 9 mm wide, in dispenser, 5.5 m or 6 m length
  - Booklet of 20 strips, 7 x 70 mm, 10 booklets per pack



Test Papers – Roll type



Test Papers – Booklet

### Ordering Information – Roll Types

Note: Type WR is 5.5 m length, others are 6 m.

Type description	Measuring pH range											Cat. No.
UNIV	1	2	3	4	5	6	7	8	9	10	11	07011030
UNIV Refill	1	2	3	4	5	6	7	8	9	10	11	07011050
WR	0	2	4	5	6	7	8	9	10	12	14	07011040
BCG	4.0	4.2	4.4	4.6	4.8	5.0	5.2	5.4	5.6		07011010	
BTB	6.2	6.4	6.6	6.8	7.0	7.2	7.4	7.6	7.8		07011020	

### Ordering Information – Booklet Types

Note: Each box contains 10 booklets of 20 strips.

Type description	Measuring pH range											Cat. No.
CR	0.4	0.6	0.8	1.0	1.2	1.4	1.6	1.8	2.0		07010010	
TB	1.4	1.6	1.8	2.0	2.2	2.4	2.6	2.8	3.0		07010020	
BPB	2.8	3.0	3.2	3.4	3.6	3.8	4.0	4.2	4.4		07010030	
PB	3.2	3.6	4.0	4.2	4.4	4.6	4.8	5.2	5.6		07010090	
PP	3.4	3.8	4.2	4.6	5.0	5.2	5.6	6.0	6.4		07010130	
BCG	4.0	4.2	4.4	4.6	4.8	5.0	5.2	5.4	5.6		07010040	
CPR	5.0	5.2	5.4	5.6	5.8	6.0	6.2	6.4	6.6		07010100	
MR	5.4	5.6	5.8	6.0	6.2	6.4	6.6	6.8	7.0		07010050	
BCP	5.6	5.8	6.0	6.2	6.4	6.6	6.8	7.0	7.2		07010140	
BTB	6.2	6.4	6.6	6.8	7.0	7.2	7.4	7.6	7.8		07010060	
PR	0.0	0.2	0.4	0.6	0.8	1.0	1.2	1.4	1.6		07010150	
AZY	10.0	10.4	10.8	11.0	11.2	11.4	11.6	11.8	12.0		07010070	
ALB	11.0	11.4	11.8	12.2	12.6	12.8	13.0	13.2	13.6		07010080	
UNIV	1	2	3	4	5	6	7	8	9	10	11	07010120

## › Litmus Paper

- **Litmus Papers** made from high purity litmus to secure high sensitivity and clear color change
- **Blue Litmus Paper** change color to red when the liquid is acidic
- **Red Litmus Paper** change color to blue when the liquid is alkaline
- **Available in two sizes and packing; booklets and loose strips**



## › Ordering Information

### Red Litmus

Size [mm]	Packing	Cat. No.
7 x 70	10 booklets each of 20 strips	07020020
9 x 85	10 boxes each of 100 strips	07020040

### Blue Litmus

Size [mm]	Packing	Cat. No.
7 x 70	10 booklets each of 20 strips	07020010
9 x 85	10 boxes each of 100 strips	07020030